
Istruzioni per l'uso / Instructions for use / Mode d'emploi / Instrucciones de uso /
Bedienungsanleitung

Europa B xp

CE
0434

I GB F E D

DXBA408 Rev_14

AUTOCLAVE MOD. EUROPA B xp

Gentile Cliente,

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AUTOCLAVE MOD. EUROPA B xp

Dear Client,

Thank you for having chosen our autoclave mod. EUROPA B xp, we know how to exchange your fidelity, with maximum attention and service definitely corresponding to your expectations.

The autoclave you have chosen is the maximum expression of applied technology. EUROPA B xp guarantees safe sterilization of all types of instruments and materials: solids, hollows, porous, free and packaged.

Before using this autoclave, we invite you to read with maximum attention the user's manual and then keep it in a place accessible to all operators in charge of sterilization.

The autoclave must be used only and exclusively for the uses foreseen by the constructor.

For installation, maintenance and assistance **ask exclusively** for a technician authorized by TECNO GAZ/MEDILINE ITALIA.

We invite you to use and ask for **exclusively** original spare parts TECNO GAZ.

TECNO GAZ SpA is a company leader in the sector and produces devices for all phases of the sterilization cycle, and can also provide **free of charge** official protocol operatives.

In accordance with the guidelines of directive 93/42/EEC

TECNO-GAZ SpA

For commercial information call : TECNO-GAZ SPA tel. 0039 0521 83 39 26

For technical information and assistance requests call: MEDILINE ITALIA srl
Tel. 0039 0522 94 29 97

Dealer where autoclave was purchased

Company name _____

Telephone _____

Technician referral name _____


TABLE OF CONTENTS

1. COMMANDS AND CONTROLS	pag.41
1.1 FRONTAL VIEW CLOSED DOOR	pag.41
1.2 FRONTAL VIEW OPEN DOOR	pag.42
1.3 BACK VIEW	pag.43
2. SAFETY	pag.44
2.1 SAFETY SYMBOLS	pag.44
2.2 SAFETY DEVICES	pag.44
2.3 SAFETY NOTICES	pag.44
2.4 DISPOSAL	pag.44
3. USE AND FUNCTION FORESEEN FOR AUTOCLAVE	pag.45
4. TECHNICAL DATA	pag.46
4.1 MECHANICAL DATA	pag.46
4.2 ELECTRICAL DATA	pag.46
4.3 CHAMBER	pag.46
4.4 CLEAN WATER TANK	pag.46
4.5 USED WATER TANK	pag.47
4.6 BACTERIOLOGICAL FILTER	pag.47
4.7 TRAY CARRIER	pag.47
4.8 TRAYS	pag.47
5. ACCESSORIES	pag.48
6. UNPACKING	pag.49
7. INSTALLATION	pag.50
8. USAGE INSTRUCTION	pag.52
8.1 TURNING ON AUTOCLAVE	pag.52
8.2 LOADING PURE WATER TANK	pag.52
8.3 MANUAL LOADING PURE WATER TANK	pag.53
8.4 CHARACTERISTICS OF WATER TO USE	pag.53
8.5 LOADING MATERIALS IN AUTOCLAVE	pag.54
8.6 BEGINNING STERILIZATION CYCLE	pag.54
8.7 END OF CYCLE	pag.54
8.8 UNLOADING STERILIZED MATERIAL	pag.54
8.9 DISCHARGING USED WATER	pag.54
8.10 INTERRUPTION OF STERILIZATION CYCLE	pag.55
9. STERILIZATION CYCLES	pag.56
9.1 CYCLES DESCRIPTION	pag.56
9.2 OPERATION CYCLES	pag.56
9.3 OPERATION CYCLES – NIGHT CYCLES	pag.56
9.4 TEST CYCLES	pag.57
9.5 CYCLE DIAGRAM	pag.57

10. MESSAGE OF ERROR OR ALARM	pag.58
10.1 LIST OF ERRORS OR ALARMS	pag.58
11. MAINTENANCE AND PREVENTION	pag.59
11.1 BACTERIOLOGICAL FILTER SUBSTITUTION	pag.59
11.2 CHAMBER CLEANING	pag.59
11.3 DOOR SEAL CLEANING	pag.59
11.4 REGULATION DOOR	pag.59
11.5 SUBSTITUTION OF PRINTER PAPER	pag.60
11.6 PRINTER STATUS	pag.61
12. PROGRAMMING AND PARAMETERS DISPLAY	pag.61
12.1 LANGUAGE SELECTION	pag.61
12.2 HOURS SELECTION	pag.61
12.3 MINUTES SELECTION	pag.61
12.4 SECONDS SELECTION	pag.62
12.5 DATE SELECTION	pag.62
12.6 MONTH SELECTION	pag.62
12.7 YEAR SELECTION	pag.62
12.8 VISUALIZE ATMOSPHERIC PRESSURE	pag.63
12.9 WATER LOADING SELECTION	pag.63
13. SOLUTIONS TO OPERATION PROBLEMS	pag.64
13.1 AUTOCLAVE DOES NOT DRY CORRECTLY	pag.64
13.2 AUTOCLAVE CHAMBER BECOMES WHITE	pag.64
13.3 AUTOCLAVE CHAMBER PRESENTS GREEN STAINS	pag.64
13.4 INTERRUPTION OF THE STERILIZATION CYCLE	pag.64
13.5 AUTOCLAVE DOES NOT RECEIVE COMMANDS	pag.64
13.6 AUTOCLAVE PRINTER DOES NOT PRINT	pag.65
13.7 STAINS ON INSTRUMENTS	pag.65
14. WATERCLEAN	pag.66
14.1 CONNECTION TO WATERCLEAN	pag.66
14.2 CONNECTION WATERCLEAN	pag.69
14.3 CHANGING WATERCLEAN RESIN CARTRIDGE	pag.69
15. STANDARD STERILIZATION PROCEDURE	pag.70
16. INSTRUCTIONS FOR AUTOCLAVE SHIPPING	pag.71
17. PROCEDURES FOR SERVICE AND ASSISTANCE	pag.72
18. GUARANTEE	pag.73
DECLARATION OF COMPLIANCE	pag. 185


1 - COMMANDS AND CONTROLS

1.1 – Frontal View Closed Door

- 1** Maximum water level
- 2** Minimum water level
- 3** Maximum used water level
- 4** Helix Test
- 5** Cycle 121°C not packaged
- 6** Cycle 121°C packaged
- 7** Cycle 134°C not packaged
- 8** Vacuum Test


- 9** Prion Cycle
- 10** Cycle 134°C packaged
- 11** Cycle 134°C Flash
- 12** Cycle Selection
- 13** Water loading
- 14** ON - OFF
- 15** Remove door lock
- 16** General switch

1.2 – Frontal View Open Door

- 1** Manual Water Loading
- 2** Bacteriological Filter
- 3** Water Loading
- 4** Used Water Disposal

1.3 – Back View

- 1 Discharge of used water
- 2 Exceeding used water – condensate discharge
- 3 Discharge clean water – connection deionizer
- 4 Exceeding clean water

- 5 Plug for Electrical power ONLY AND EXCLUSIVELY FOR DEIONIZER
- 6 Safety valve
- 7 Fuse 12A
- 8 Fuse 12A
- 9 Fuse 500mA
- 10 Connection RS232
- 11 Principal electrical alimentation

2 - SAFETY

2.1 – Safety Symbols

HIGH VOLTAGE

SHUT OFF VOLTAGE BEFORE REMOVING COVER

HIGH TEMPERATURE

CONNECTION TO GROUND

2.2 Safety Devices

Safety devices are the following:

- Safety valve calibrated 2.4 bar $\pm 10\%$

- Electromagnetic lock to avoid opening of door during cycle

2.3 Safety Notices

In quality of constructor Mediline Italia Srl is responsible of the product commercialised, to senses of the enforced laws . The responsibility expires at the moment the device, or parts of the device, is tampered by a non authorized staff or if the parts are not original.

The location where the autoclave is installed must not be at risk for potential explosion and/or fire.

The autoclaves must be installed in an atmosphere in compliance with legislative requirement in force on the Country use.

2.4 Disposal

Disposal of the autoclave must be effected according to current legislation in the country of use. No responsibility is assumed for disposal operations executed in contrast to current legislation.

3 - USE AND FUNCTION FORESEEN FOR AUTOCLAVE

Autoclave mod. EUROPA B xp is able to sterilize the 3 kind of charge indicates in the European Norm EN 13060, specifically :

FERROUS MATERIALS OR SOLIDS Instruments or devices without cavity or obstacles for steam penetration

POROUS PARTS Simple or composed materials they can absorb the fluid (textures, overalls, gauzes, medications, etc...)

HOLLOW PARTS Materials or devices with cavity, obstructions, etc...
European Norm divides these kind of parts in 2 types, it depends by dimensions of cavity. Indicatively

TYPE B : tubes with large hollow

TYPE A : turbines, devices and strick with small holes or small dimensions

Autoclave must be used, only and exclusively, for the sterilization of instruments and materials compatible with steam sterilization systems.

This equipment can be used in the dental, medical, aesthetic fields and, generally, in all the fields where the sterilization of the instruments and materials is made

This device has been certified for the sterilization of the following materials:

-wrapped and unwrapped solid materials for a maximum load of 5 Kg.

-Type A/B hollow materials for a maximum load of 5 Kg

-Porous materials for a maximum load of Kg.1,5

Valid for European countries only

Always verify that the load to be sterilized can support the cycle temperature chosen.

4 – Technical Data

4.1 – Mechanical Data

- Working temperature +5 +40°C
- Maximum altitude 2.000 m
- MAX relative humidity at 30°C 80%
- MAX relative humidity at 40°C 50%
- Dimensions of space occupied (L x H x P) 510x390x590 mm
- Space occupied with open door 370 mm
- Weight (tank empty) 54 Kg
- Weight (tank full) 63 Kg
- Weight of area of support 2058 N/m²
- Volume 0.12 m³
- Color RAL 9016
- Material Stainless steel AISI 304/ Steel FeP01
- Potential sound level < 70 db A

4.2 – Electrical Data

- Power voltage 230 V a.c. + \- 10 % single phase
- Power 1700 W
- Frequency 50 / 60 Hz
- Power cord 2 + 1 x 1mm²
- Fuses 6.3 x 32 - 12 A
- Heat transmittend 3.6 E +6 J \ hour

4.3 - Chamber

- MAX working pressure 2.2 bar (relative)
- MAX empty - 0.9 bar (relative)
- MAX Temperature 138 °C
- Material Stainless steel AISI 304
- Dimensions diam. 245 x 318 mm

4.4 – Clean Water tank

- Volume 4.5 l
- Usable cycles 5
- Material Polyethylene

4.5 – Used Water tank

- Volume 4.5 l
- Usable cycles 5
- Material Polyethylene
- Max temperature used water 50°C

4.6 – Bacteriological filter

- Diameter 56 mm
- Filtering capacity 0.3 micron

4.7 – Tray carrier

- Material Anodized aluminum
- Dimensions (L x H x P) 192 x 165 x 280 mm
- Number provided 1

4.7.1 – Tray carrier optional

- Material Anodized aluminum
- Dimensions (L x H x P) 193 X 200 X 280) mm
(200 X 193 X 280) mm

4.8 - Trays

- Material Anodized aluminum
- Dimensions (L x H x P) 184 x 17 x 286 mm
- Number provided 4

5 - ACCESSORIES

Autoclave is sold and delivered complete of all accessories, which are included in the original packaging.
Accessories provided:

- N. 1 PC. CLOSED TRAY CARRIER IN ANODIZED ALUMINUM
- N. 4 PC. PERFORATED TRAYS IN ANODIZED ALUMINUM
- N. 1 PC. KEY FOR EXTRACTION TRAY AND REGULATION DOOR
- N. 1 PC. PLASTIC FUNNEL WITH TUBE
- N. 1 PC. HOSE FOR LOADING WATER WITH THE FILTER
- N. 1 PC. HOSE FOR DISCHARGING WATER
- N. 1 PC. RILSAN HOSE
- N. 1 PC. SPONGE FOR CLEANING CHAMBER

ATTENTION: NO WARNINGS ARE NEEDED FOR THE ACCESSORIES ABOVE

IMPORTANT

Ask for and use only and exclusively original accessories Tecno - Gaz.

6 - UNPACKING

The autoclave is shipped in a suitable package to be transported and moved easily and to protect its contents. The package must not be jostled, it must be moved with care avoiding letting it turn over or fall from a height superior to 15 cm.

In case autonomous handling means are not available handle the packaged autoclave always in two persons.

The autoclave is supported on a wood pallet and closed in a corrugated reinforced carton, inserted from above by machine.

The autoclave is fixed with two plastic bands to the pallet.

To unpack the autoclave, cut the two bands, then remove the carton covering it with a vertical movement.

The autoclave will be resting on the pallet, lift it always in two persons, picking it up at the base and positioning it on the prepared surface or in the room.

IMPORTANT

Do not lift the autoclave from the inferior part of the door or command panel, this incorrect operation can create problems of a mechanical nature.

Inside the package you will find:

USER'S MANUAL to be read with attention and kept in a place available to all operators assigned to sterilization.

CERTIFICATION which must be conserved

GUARANTEE LABEL which must be completed and sent to TECNO GAZ SPA.

ATTENTION

ALWAYS CONSERVE ORIGINAL PACKAGING


7- INSTALLATION

- 1 – Install the autoclave in a suitable ambient for sterilization.
- 2 – The room must be adequately lighted and aired, as foreseen in the current directives.
- 3 – Install the autoclave far from heating sources and sprinkles of water.
- 4 – Position the autoclave on a suitable surface to support the weight and of adequate dimensions.
- 5 – The supporting surface must be perfectly horizontal and must not have inclinations.
- 6 – Position the autoclave at a height where the user can inspect the entire sterilization chamber and clean it easily.
- 7 – Open the door of the autoclave and remove from inside the sterilization chamber all packages which contain the single accessories.
- 8 – Leave inside the sterilization chamber only the tray carrier with the trays. All other accessories should be positioned in a separate space available to operators.
- 9 – Never place on the autoclave newspapers, trays, liquid containers, etc.
- 10 – Do not lean on the door when it is open.
- 11 – Please keep at least 15 cm of free space from the back and from the side of the autoclave to achieve the correct ventilation.
- 12 – Always check that the electrical system to which the autoclave is attached is in compliance with current regulations and dimensioned in a an adequate manner for the apparatus's characteristics.
- 13 – Take the electrical power cord provided and insert the female plug in the socket of the back panel of the autoclave (fig.1 pos.6).
Insert the electric plug of the system checking that the power is 230 Vac- 2000Va.

IMPORTANT

Avoid connections with extension cords, reductions or adapters, this can create micro interruptions with consequent alarm indications.

FIGURE 1


Back view 1.3
Pos. 6

14 - Turn on the autoclave by pushing the general switch (fig.2 pos.16) and open the door to the autoclave. Wait a few seconds, there will be two beeps which indicate the reaching of relative parameters for automatic atmospheric alignment, at the same time the display will indicate "door open".

IMPORTANT

Never select a command before the two beeps, in that the autoclave will not accept the chosen program.

FIGURE 2

Frontal view 1.1
Pos. 16

15 - Load the pure water tank following the indications in the section entitled USAGE INSTRUCTIONS. Successively run the test cycle, always following the indications expressed in the section USAGE INSTRUCTIONS.


16 - If the autoclave is installed with a system of deionization mod. WATERCLEAN, comply with instructions contained in the package of this device.

8- USAGE INSTRUCTION

After having installed the autoclave, proceed with preparations and use.

8.1 – TURN ON THE AUTOCLAVE, press the general switch (fig.3 pos.16)

FIGURE 3

Frontal View 1.1
Pos. 16


8.2 – LOADING PURE WATER TANK (lt. 4.5)

Connect the transparent hose complete with the filter, (provided), to the frontal connection of the autoclave (fig.4 pos.3); unscrew the metal ring anticlockwise (to execute 2 complete turns)
Insert the other end of the hose with filter to the inside of the demineralized or distilled water container.

At this point press the button  to start the water loading pump and press it for at least one second.
The pump loads the tank inside the autoclave. If maximum level is not reached within 180 seconds the pump automatically stops, therefore it is necessary to press the button again.
The pump stops automatically when maximum level is reached.

With connection to a deionizer, if the maximum water level has not been reached, the functioning of the autoclave will be obstructed.

FIGURE 4


Frontal view 1.2
Pos. 3

8.3 – MANUAL LOADING OF PURE WATER TANK

In the case the water loading pump does not function, filling the tank can be done manually by the operator using the following method:


- turn on autoclave
- twist off the top located on the front of the autoclave (fig. 5 pos.1)
- insert in this opening the hose connected to the funnel
- pour distilled water in the funnel, keeping it higher than the loading top
- add until the led of maximum level turns on (fig.6 pos.1)
- loading completed, twist on the top

FIGURE 5

Frontal view 1.2
pos. 1

FIGURE 6

Frontal view 1.1
pos. 1

8.4 – CHARACTERISTICS OF WATER TO USE

Attachments

TABLE OF QUALITATIVE LEVELS ESTABLISHED BY REGULATION DIN EN 285

CEN STANDARD DIN EN 285			
Evaporation residual	≤	10	mg/l
Silicon oxidation	≤	1	mg/l
Iron	≤	0.2	mg/l
Cadmium	≤	0.005	mg/l
Lead	≤	0.05	mg/l
Residual of heavy metals such as iron, cadmium, lead	≤	0.1	mg/l
Chloride (Cl ⁻)	≤	2	mg/l
Phosphate (P20s)	≤	0.5	mg/l
Conductivity (at 20°C)	≤	15	μs/cm
PH value (grade of acidity)	5 to 7		
Appearance	Colorless, limpid, without deposits		
Hardness (E alkaline-earth ions)	≤	0.02	mmol/l

8.5- LOADING MATERIAL IN AUTOCLAVE

Place the materials to be sterilized on trays provided.

-do not pile on each other the materials

-place instruments in envelopes always with the paper part towards the top

-never place the materials in direct contact with the sterilization chamber or the closed door

Upon finishing loading, close the door of the autoclave. On the display DOOR CLOSED will appear.

8.6 – BEGINNING OF STERILIZATION CYCLE

Choose sterilization program most adaptable of load prepared, press the button

Each individual pressure changes program selection to the next.

After having chosen the program start the cycle by pressing the button

During the cycle, the display will show all parameters and relative information.

8.7 – END OF CYCLE

A beep will indicate to operators the finished sterilization cycle.

If the cycle has been successful, the display will show the message END CYCLE, in case of problems, it will show a code of error or alarm.

When the acoustic signal finishes, unlock the door pressing the button

Inside the chamber, the button will not function the unlocking. Wait until the chamber is completely depressurized and repeat the operation or, only in cases of emergency, press simultaneously Door + Pump. When the door is unlocked, turn the door handle and open.

8.8 – UNLOADING STERILIZED MATERIALS

Wear protective gloves, as requested by that is obligatory for the operator security.

Extract the tray using the suitable key provided, place the instruments or full trays in closed and protected spaces.

8.9 – DISCHARGING USED WATER

When the led of used water level (fig.7 pos.3) goes on, proceed with emptying. **If not emptied, the autoclave function is inhibited.**

Take the hose provided and insert it in the used water discharge plug, on the front of the autoclave (fig.8 pos.4), unscrew the metal ring anticlockwise (to execute 2 complete turns), simultaneously place the other end of the hose in a container, the water will descend into the same container.


IMPORTANT

a – the hose located in the receiving container, must not ever, touch or be immersed in the discharged water, otherwise there will be a situation of sucking up.

b- always wait until the unloading water is totally emptied. The led of maximum used water level will turn off when there is still water in the tank, therefore do not use as a reference for this operation.


To the term of the drainage to screw the metal ring and to remove the tube.

FIGURE 7

Frontal view 1.1
Pos. 3

FIGURE 8


Frontal view 1.2
Pos. 4

8.10 – INTERRUPTION OF STERILIZATION CYCLE

A sterilization cycle can be voluntarily interrupted by pressing the button  for at least 2 seconds. The autoclave will emit a beep, go into decompression and on the display will appear the message error 001 (INTERRUPTED CYCLE).

To reset error, maintain pressed simultaneously the buttons below indicated, until the intervention of an acoustic signal:

9 – STERILIZATION CYCLE

9.1 – CYCLE DESCRIPTIONS

The autoclave has three series of cycles:

- a- operation cycles
- b- operation cycles-night cycles
- c- test cycles

9.2 – OPERATION CYCLES

All operation cycles have a system of fractioned vacuum, which therefore permits sterilization of materials that are hollow, porous, solid, both free and packaged.

Temperatures can be selected from 121°C – 134°C.

Normally the cycles of 121°C are used for thermoplastics or sensitive materials, while the 134°C cycles are used for all other materials.

In all cases always follow the indications given by the manufacturer of the instruments or devices to be sterilized.

There are also two special cycles:

-prion cycle, studied for the Crutzfeldt-Jacobs disease (mad cow disease)

-flash cycle created for rapid sterilization of unpackaged instruments and devices.

9.3 – OPERATION CYCLES – NIGHT CYCLES

Europa B xp is provided with a special economizer device.

This device allows the execution of sterilization cycles without the presence of an operator.

At cycle end, if the door is not opened, the autoclave stabilizes and then shuts off automatically. Only the general switch remains on.

Upon the operator's arrival it is sufficient to press any button to re-start the autoclave and read the cycle result on the display.

Besides that the printer will have regularly prepared the written report relative to the same cycle.

STERILIZATION TABLE

CYCLES	TIME EXPOSED (Minutes)	DRYING TIME (Minutes)	CHARGE* (Kg)	CHARGE** (Kg)	MINIMUM WORK PRESSURE (bar, differential pressure)	VACUUM PHASES N°
121° Unpackaged	18	6	1.25	0.375	1.04	3
121°C Packaged	18	12	1.25	0.375	1.04	3
134°C Flash	4	2	1.25	0.375	2.02	3
134°C Unpackaged	5	6	1.25	0.375	2.02	3
134°C Packaged	5	12	1.25	0.375	2.02	3
134°C Prion	20	12	1.25	0.375	2.02	3

* The values indicate are the maximum load for every tray, for ferrous materials or solid or hollow part

** The values indicate are the maximum load for every tray, for porous parts

Is not possible to add the maximum loads indicates in the table over.

Warmig up time and vacuum fractionation can vary, because of autoclave's conditions, between 25 and 30 minutes. You have to an this time to sterilization's times showed in table

9.4 – TEST CYCLES

-Helix test: this test serves to verify steam penetration capacity in a hollow body. Use the Helix Indicator System.

- Bowie and dick test: this serves to verify steam penetration capacity in a porous load. Use Bowie and Dick test pack (top table).


- Vacuum test: this serves to verify the preservation capacity of the autoclave vacuum; test report is prepared by the autoclave printer.

It is preferable to use tests distributed by TECNO GAZ spa.

TEST CYCLES TABLE

CYCLES	STERILIZATION TIME (Minutes)	DRYING TIME (Minutes)	MAX. LOADING FOR TRAY (Kg)	MINIMUM WORK PRESSURE (bar, differential pressure)
134°C Helix/ Bowie&Dick	3.5	4		2.02
134°C Vacuum				-0.9


CYCLE DIAGRAM

10 – MESSAGES OF ERROR OR ALARMS

Messages of error are emphasized through an Alfa-numerical code, consisting in a letter and three numbers. The “E” suffix is relative to errors of the operator, which can be corrected by same operator. The “A” suffix is relative to alarms, irregularities of the autoclave, after having carried out the remedy, if the problem persists it is necessary to request telephonic assistance.

To reset alarms and errors, maintain pressed the buttons here indicated until the intervention of a beep:

10.1 – LIST OF ERROR-ALARM MESSAGES

ERROR	CAUSE	REMEDY
E 200	Low water level	Fill tank
E 250	Low level water cleaned up with loading from deionizer	Fill tank
E 300	Maximum level used water	Empty tank
E 401	Door open	Close door with attention
E 451	Door unlocked	Open door and close again
E 461	High temperature in chamber for Vacuum test (above 40°)	Shut off and wait 10 minutes
A 250	Cycle begins without water	Control Waterclean device
A 001	Interrupted cycle	Reset and restart cycle
A 011	Irregular pc display	Shut off and restart autoclave
A 101	Vacuum not reached in 10min.	Reset and repeat the cycle
A 111	Vacuum not maintained for first phase VT	Reset and repeat the cycle
A 121	Vacuum not maintained for second phase VT	Reset and repeat the cycle
A 401	Hatch opened during the cycle or problem closing	To control the closing system
A 501	Pressure not maintained	Reset and repeat the cycle
A 551	Irregular 1° temperature sensor	Reset and repeat the cycle
A 552	Irregular 2° temperature sensor	Reset and repeat the cycle
A 641	Overly high pressure	Reset and repeat the cycle
A 651	Overly high temperature	Reset and repeat the cycle
A 661	Irregular temperature comparison	Reset and repeat the cycle
A 701	Working pressare not reached	Reset and repeat the cycle
A 751	Low temperature	Reset and repeat the cycle
A 761	Irregular pressure sensor	Reset and repeat the cycle
A 801	Irregular unloading pressure	Reset and repeat the cycle
A 901	Main voltage too low	Inform electric energy provider
SUBSTUTUTE RESIN CLEN WATER resin finished		Follow the instruction in Paragraph 14.2 of the manual
D	The printer's door is not closed	Try to close the printer's door
P	Paper is missing	Insert paper
L	There is low volatge on the line	To try to connect the machine to the net electrical worker in an other point; if the problem persists to contact a electrician

IMPORTANT

If after following the remedies, the problem persists, call telephone assistance at MEDLINE ITALIA srl 0039 – 0522 – 94 29 97.

11 – MAINTENANCE AND PREVENTION


Correct maintenance of the autoclave assures correct functioning of it and a secure saving in terms of time and costs for assistance and maintenance.

11.1 – BACTERIOLOGICAL FILTER SUBSTITUTION

Substitute the bacteriological filter (Fig.9 pos.2), maximum every 200 cycles, or as soon as the filter becomes a dark color.

Use only and exclusively original TECNO GAZ filters.

FIGURE 9

Frontal view 1.2

Pos. 2

11.2 – CHAMBER CLEANING

Clean periodically the chamber, remove eventual deposits or debris, thus avoiding the introduction in the discharging circuit of material which can cause obstructions.

For correct cleaning **use only water** and the abrasive sponge provided.

Never use solvents, detergents, chemical solutions, descaling agents or other similar products.

11.3 – DOOR SEAL CLEANING

Periodically remove eventual residuals which are deposited on the circumference of the seal using water and the sponge provided (non abrasive side), or else a moistened cloth.

11.4 – DOOR REGULATION

To regulate the door lock pressure, it is necessary to adjust the regulator of the same door using the double function wrench (extraction tray + wrench) provided.

Turn anti clockwise (fig.A) to increase lock pressure.

Turn clockwise (fig.B) to decrease the lock pressure.

FIG. A




FIG. B

11.5 – SUBSTITUTION PRINTER PAPER

Thermographic paper roll 57 mm wide. Thermographic paper can be printed on one side only. Open the printer cover (1), place the paper roll by paying attention to the rotation direction (2), pull the paper until it comes out and close the cover (3); the printer is ready for printing (4).

11.6 – PRINTER STATUS

If the cover of the printer is not perfectly closed, on the display the letter "D" appears.
Open and close the printer's cover.

If paper is missing, on the display the letter "P" will appear.
Open the cover, insert paper and re-close the cover.

12 – PROGRAMMING AND PARAMETERS DISPLAY

12.1 – Language Selection

Press simultaneously  and MAINS switch. Hold down 

until on the display appears the name of the language programmed.

Press  to modify the language (forward).

Press  to modify the language (backwards)

12.2 – Hours Selection

Press  the hours appears.

Press  to modify hour (forward).

Press  to modify hour (backwards).

12.3 – Minutes Selection

Press  minutes appear.

Press  to modify minutes (forward).

Press  to modify minutes (backwards).

12.4 – Seconds Selection

Press  seconds appear. Do not modify.

12.5 – Date Selection

Press  day appears (flashing).

Press  to modify day (forward).

Press  to modify day (backwards).

12.6 – Month Selection

Press  month appears (flashing).

Press  to modify month (forward).

Press  to modify day (backwards).

12.7 – Year Selection

Press  year appears.

Press  to modify year (forward).

Press  to modify year (backwards).

12.8 - Visualize Atmospheric Pressure

Press  memorized atmospheric pressure appears. **DO NOT MODIFY!**

12.9 – Water Loading Selection

Press  water loading position appears (default: **loading with pump**).

Press  to modify to water **supply loading**

Press  and maintain pressed until the sound of the “beep”.

Press  and maintain pressed until the sound of the “beep”.

The data inserted are memorized. Turn off the autoclave using the MAINS switch.

13 – SOLUTIONS TO OPERATIONAL PROBLEMS

In many cases some alarms or errors are determined by lack of attention or lack of familiarity with some technical and operational aspects.

Following are listed some cases of irregularity with relative solutions.

13.1 – AUTOCLAVE DOES NOT DRY CORRECTLY

a - SUBSTITUTE THE BACTERIOLOGICAL FILTER, with a new TECNO GAZ original.

b - ORIGINAL TRAYS HAVE NOT BEEN USED, but other trays without perforations or diverse perforations. Use only original trays, possibly in aluminum. Request additional original trays.

c - INSTRUMENTS HAVE NOT BEEN ARRANGED CORRECTLY, arrange instruments or packaged materials with the paper turned upwards.

Never pile up instruments or materials, arrange them parallel to one another.

Never place the instruments or materials in contact with the chamber, they can burn or become deformed.

13.2 – AUTOCLAVE CHAMBER BECOMES WHITE

a – CHANGE IMMEDIATELY THE TYPE OF WATER USED, use demineralized or distilled water, as is specifically indicated in the previous chapter and then proceed with cleaning the chamber.

b – THE WHITISH COLOR CAN be a consequence of the evaporation of organic materials present on the instruments.

Instruments should be cleaned with a detergent action more suitable and thoroughly.

13.3 – AUTOCLAVE CHAMBER PRESENTS BLUEISH GREEN STAINS

a – INSTRUMENTS HAVE NOT BEEN CORRECTLY RINSED cleanse and rinse with major attention the instruments.

If the stains are evident request telephonic assistance.

13.4 – STERILIZATION CYCLE INTERRUPTS WITHOUT APPARENT MOTIVE

a – CONTROL IF AUTOCLAVE IS CONNECTED to the electrical net with extensions, reductions, adapters, if so remove this accessories and connect the autoclave directly to the electric socket.

13.5 – AUTOCLAVE DOES NOT RECEIVE COMMANDS

a – AUTOCLAVE PRE-ARRANGES BAR AUTOMATIC ALIGNMENT wait for the double beep after opening the door, then program the functions.

b – PURE WATER TANK is empty, minimum led level is on, fill with pure water

c – USED WATER TANK IS FULL. Maximum level led is on, discharge used water.

13.6 AUTOCLAVE PRINTER DOES NOT PRINT

- a** - PAPER CARRIER COVER IS NOT CLOSED CORRECTLY, on the display appears the letter "D". Open and close the printer cover checking that is correctly closed.
- b** - PAPER MISSING, on the display appears the letter "P". Open the cover and insert the roll of thermo-graphic paper.
- c** - ROLL OF PAPER IS MOUNTED WITH THERMO SIDE TOWARDS INSIDE, open the cover and turn the roll of paper to the correct position.

13.7 STAINS ON INSTRUMENTS

- a** - INSTRUMENTS BECOME YELLOW, chemical liquid residual has with heat been fixed on the instruments. They were not sufficiently rinsed.
- b** - STERILIZATION CHAMBER PRESENTS YELLOW SPOTS, instruments were introduced in the chamber with chemical liquid which dripped and is fixed due to heat. Rinsing was not sufficient.
- c** - INSTRUMENTS PRESENT WHITE SPOTS, rinsing was done with water with high level of calcium and the instruments were not dried. For last rinse it is advisable to used demineralized water and dry carefully the instruments.
- d** - INSTRUMENTS ARE DARKENED, this is caused by the fact that the instruments have internally a strong component of carbon.

14- WATER CLEAN

14.1- CONNECTION TO WATER CLEAN

Europa b xp autoclave is pre-arranged for loading demineralized water through the demineralizer apparatus WATER CLEAN.


WATER CLEAN is a demineralizer apparatus , accessory to TECNO GAZ S.p.A. autoclave (not included). This product is used to eliminate mineral salts in water coming from the water supply to feed the autoclave. Before the operator installs WATER CLEAN he must program the autoclave following the instructions here below:

1.	Switch-on the autoclave with the main (green) switch and keep pressed  (during the swith-on) until you see a string reporting the selected language.
2.	Press and relase  to select next screen until you see: - charging by pump or - charging by water conditioner To change this settino use 
3.	To save keep pressed  until you hear a "BEEP"
4.	To exit keep pressed  until you hear a "BEEP"

The diagram illustrates a water treatment system for an autoclave. The system components and their connections are as follows:

- AUTOClave**: The main unit being treated, featuring a control panel with a digital display and buttons (1, 2, 3, 4, 5, 6, 7, 8, 10, 11).
- Water Clean Unit**: A unit that treats water, connected to the autoclave via a line with valve **SWCA045**.
- Water Supply**: A line with valve **SWCA019** that provides water to the Water Clean Unit.
- Control Valves and Sensors**:
 - SWCA046**: A valve located on the line between the Water Clean Unit and the autoclave.
 - SWCA022**: A valve located on the line between the Water Clean Unit and the autoclave.
 - SWCA045**: A valve located on the line between the Water Clean Unit and the autoclave.
 - SWCA047**: A valve located on the return line from the autoclave.
 - Sensors**: 12A, 12B, 12C, 43, and 38 are located on the water supply and treatment lines.
 - Sensor 24**: A sensor located on the return line from the autoclave.

SCHEME OF CONNECTION OF THE VALID DEMINERALIZER FROM THE 12/05

14.2 – CONNECTION OF WATER CLEAN

Turn off the autoclave by pressing the General switch located on the right side of the autoclave. The operator after having set the autoclave can install the demineralizer apparatus WATER CLEAN following the indications in the manual attached to the demineralizer apparatus

Here are indicated the specific connection to the autoclave of the water supply hose and the connection of the electric plug of command of WATER CLEAN:

5. Turn off the autoclave if it is on;
6. Close the faucet located above the demineralization system;
7. Install WATER CLEAN as indicated in its manual;
8. Twist the male screw of the connection-hose carrier SWCA047 with Teflon or other component which guarantees the sealing of water;
9. Screw the connection-hose carrier SWCA047 on the female screw to the unloading clean water pos.(3);
10. Insert the SWCA045 hose exit of CLEAN WATER to the connection-hose carrier SWCA047 already screwed to the autoclave
11. Seal with a band the hose SWCA045 to the connection-hose carrier SWCA047 with a band to avoid water loss;
12. Insert the plug Pos. (24) of WATER CLEAN to the socket POS. (5) on the back of the autoclave
13. Open the faucet located above the demineralization system;
14. Control that there is no water loss;
15. Turn on the autoclave;
16. Run one or more sterilization cycles to control the functioning of the connection and check above all for leaking.

AT THE END OF EACH DAY ALWAYS CLOSE THE FAUCET LOCATED ABOVE THE DEMINERALIZATION SYSTEM TO PREVENT FLOODING

CONNECT WATER CLEAN ONLY TO TECNO-GAZ PREPROGRAMMED AUTOCLAVE .

IF THE UNIT AT THE BEGINNING OF THE CYCLE SHOWS E250 THEN YOU HAVE TO FILL THE CLEAN WATER TANK UNTIL THE MAXIMUM IS REACHED. (TO WATCH PAR.10)

NOTE: FOR THE CONNECTION OF WATER CLEAN TO THE AUTOCLAVE TECNO-GAZ; REFER TO INDICATIONS IN THE WATER CLEAN MANUAL.

14.3 – CHANGING WATER CLEAN RESIN CARTRIDGE

Autoclave display indicates the state of finishing of the resin with the message "SUBSTITUTE RESIN". Order from a TECNO-GAZ supplier the resin cartridge code 242-S; Substitute the resin cartridge as indicated in the WATER CLEAN manual.

After substitution, press button  until the "beep", to reset the autoclave.

15 – STANDARD STERILIZATION PROCEDURE

Sterilization means adopting precise operational protocol and applying a system.

Autoclave is the final link of this system.

Those who do not follow all of the various phases of the sterilization system will be unsatisfied with the final result.

TECNO-GAZ SpA produces all the devices for each phase of the sterilization cycle.

Therefore it is necessary to:

1 - Prepare a precise written operational protocol to distribute to all operators assigned to sterilization.

TECNO-GAZ can provide these personalized protocols.

2 - Follow carefully and with attention all of the sterilization cycle phases:

DISINFECTION obligatory phase, in base of that is obligatory for the operator security, to be done with immersion in chemical liquids or thermo-disinfection.

Ref. TECNO GAZ - vasche mod. artica – termodisinfettore mod. aqua

CLEANSING the most important phase which insures the removal of all types of residual, chemical and organic. The most suitable instruments are ultrasound baths.

Ref. TECNO GAZ – ultrasound mod. tecna

DRYING essential phase to avoid corrosion of instruments and interference of the sterilization cycle.

Ref. TECNO GAZ – device mod .unika

PACKAGING essential phase for maintaining sterile the treated instruments.

Ref. TECNO GAZ – sealer mod .opera

STERILIZATION final phase of steam sterilization

Ref. TECNO GAZ –autoclave mod. europa b xp


16 – INSTRUCTION FOR AUTOCLAVE SHIPPING

In case the autoclave must be shipped for repairs, controls, reactivation, revisions, validations follow the obligatory indications below:

- 1 – Use the original packaging; if this is no longer in your possession, use adequate packaging. The merchandise travels at risk to the sender.
- 2 – Ship only the autoclave, always removing trays and tray carrier.
- 3 – Carefully clean the sterilization chamber and autoclave in general before shipping. In case it arrives dirty or with residual the autoclave will be returned without being repaired, or it will be put through a cleansing action and disinfection.
- 4 – Always empty the clean water tank through the attachment located on the back of the autoclave (fig. 10 pos.3).
- 5 – Always empty the used water tank through the attachment located on the back of the autoclave (fig. 10 pos. 4).

FIGURE 10

Frontal view 1.3
Pos. 3

6 – Indicate by letter and insert in the package a document which indicates precisely the irregularity or service desired.

7 - Ship at your expense, otherwise you will be billed for shipping.

ATTENTION

All non original packaging which arrives will be disposed of.

Autoclave will be returned with new and original packing to insure maximum protection for your autoclave during shipping.

Cost of packaging will be charged to client.

17 - PROCEDURE FOR SERVICE AND ASSISTANCE

Repairs – revisions – assistance – validation

We consider it important to inform you of the procedures used by our service department.

The service department is the autonomous entity MEDILINE ITALIA SRL (TECNO-GAZ GROUP), which has an articulated and professional structure.

Procedures applied have been explained in our quality manual.

1 - When the autoclave arrives at MEDILINE ITALIA SRL it is consigned directly to the service center, which codes and creates the technical assistance document.

2 – Within 48 hours, it is examined by technicians and technical analysis along with cost estimate is prepared for the service requested.

If the autoclave is not delivered as requested in the chapter INSTRUCTIONS FOR SHIPPING, above times can sensibly differ

3 – Cost estimate will be sent to the distributor client, who will transmit it to the final client for vision and approval.

4 – After MEDILINE ITALIA SRL receives written acceptance of the cost estimate, the autoclave will be serviced and reshipped according to times indicated on the cost estimate.

5 – The merchandise will be reshipped in original packaging by the client's shipper.

18 – GUARANTEE

1. The product is guaranteed for a period of five years excluding electrical and electronics parts for which the guarantee is 12 months.
2. The guarantee begins from the date of delivery of the machine to client proven by the restitution of the guarantee label attached to the instructions for installation, filled in, stamped and signed by the dealer.
3. In case of contestation, the date indicated on the sales receipt is considered valid.
4. Repairs or substitutions during guarantee of a part is effected upon the absolute judgment of TECNO-GAZ company, and does not include the transfer of employees, packaging and shipping costs.
5. Excluded from the guarantee are light bulbs, fuses, and breakdowns and damage caused by improper maintenance, negligence, incompetence or other causes not related to the constructor.
6. Components subject to normal wear and accessories (for ex. buttons, parts in movement of pumps, etc) are not included in the guarantee.
7. The right to complete substitution of the machine is not recognized.
8. The guarantee does not include any compensation for direct or indirect damage of any nature towards persons or things, caused by the inefficiency of the machine.
9. TECNO-GAZ does not respond to breakdowns or damage caused by inappropriate use of the product or on that where ordinary maintenance has not been practiced or that elementary principals of correct maintenance (negligence) have been ignored.
10. There is no acknowledgment of compensation for machine stoppage.
11. The guarantee expires automatically when the machine is tampered with, repaired or modified by the buyer or by third parties not authorized by TECNO-GAZ.
12. For reparations, the buyer must consult only the dealer or the service centers indicated by TECNO-GAZ.
13. Components substituted in guarantee must be returned to TECNO-GAZ postage paid.
14. Without restitution the cost of part will be charged to the buyer.
15. TECNO-GAZ does not accept return of final users.
16. Return to TECNO-GAZ for repairs, must be controlled by the dealer or by service centers chosen by the final user in compliance with commercial procedures CM-P-003 entitled "CONTROL OF REQUESTS FOR RETURN, REPAIR, SUBSTITUTION OF TECNO-GAZ PRODUCTS AND CLIENT COMPLAINTS".
17. Return to TECNO-GAZ must be documented and authorized by internal procedures.
18. Products returned to TECNO-GAZ must have attached documentation of authorized return and a document describing the breakdown.
19. All products in repair must be shipped to TECNO-GAZ postage paid and adequately packaged (we suggest using the original packaging).

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TECNO-GAZ S.p.A.